

## SG NON-INDEXED PCR-FREE LIBRARY PREP ADAPTER OVERVIEW

The SG Non-Indexed PCR-Free Library Prep Adapter kit provides adapter and reagents designed to support the preparation of non-indexed PCR-free libraries for the G4™ Sequencing Platform.

For more information about Singular Genomics® primers and adapters, see the [Adapters and Indices for the G4 Sequencing Platform Reference Guide](#).

### Product Configuration

Product Name	Ordering #	Quantity
SG Non-Indexed PCR-Free Library Prep Adapter (24 Rxns)	#700,115	24 reactions

### SG Non-Indexed PCR-Free Library Prep Adapter Contents

One SG Non-Indexed PCR-Free Library Prep Adapter product contains:

Components	Container	Description
SG Non Indexed PCR-Free Library Prep Adapter	Freestanding micro-centrifuge tube	Adapter for PCR-free library preparation containing S1, SP1, S2, and SP2 sequences without indices, 15 µM in 130 µL.
4X Cleave Reaction Buffer	Freestanding micro-centrifuge tube	Buffer for cleave reaction. 920 µL per tube.
Cleave Enzyme	Freestanding micro-centrifuge tube	Enzyme that cleaves the cleavable site in the adapter. 24 rxn: 130 µL.

### Consumable Storage Temperatures

Components	Shipping Temperature	Storage Temperature
SG Non-Indexed PCR-Free Library Prep Adapter	Dry ice	-15°C to -25°C
4X Cleave Reaction Buffer	Dry ice	-15°C to -25°C
Cleave Enzyme	Dry ice	-15°C to -25°C

### Shelf Life

Always check the expiration date on your kit before use. Expiration dates are listed on the kit label.

For Research Use Only. Not for use in diagnostic procedures.

## Prepare SG Non-Indexed PCR-Free Library Prep Adapter for Use

1. Remove adapter and 4X Cleave Buffer from -25°C to -15°C storage.
2. Place the adapter and 4X Cleave Buffer in a water bath with room temperature deionized water.
3. Leave reagents until completely thawed.
4. Vortex thawed reagent tubes and then centrifuge briefly.
5. Keep on ice until ready to use.

**NOTE**

Do not thaw or vortex the tube with Cleave Enzyme. Keep in -25°C to -15°C storage as long as possible, and when taken out of storage, keep on ice or in a -25°C to -15°C cooling block.

### Revision History

Document #	Revision	Release Date	Description of Change
600018	Rev. 0	June 2023	Initial release.

## Customer Care

Our team is committed to ensuring your success. Field application scientists, field service engineers, and customer care specialists are dedicated to responding quickly and effectively to your inquiries. From preparing your lab for the arrival of your instrument, to onsite training in system operation and maintenance and continued support, we'll be with you every step of the way.

<b>Website</b>	<a href="https://singulargenomics.com">singulargenomics.com</a>
<b>Email</b>	<a href="mailto:care@singulargenomics.com">care@singulargenomics.com</a>
<b>Telephone</b>	442.SG.CARES(442.742.2737)

Singular Genomics Systems, Inc.

3010 Science Park Rd

San Diego, CA 92121 U.S.A.

This product is covered by one or more patents, trademarks and/or copyrights owned or controlled by Singular Genomics Systems, Inc.

Singular Genomics and G4 are trademarks of Singular Genomics Systems, Inc.